AQUA

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RESULTADOS	TEST DE PROCEDIMIENTO															
<p>Bacterial Count</p>  <p>Low contamination Moderate contamination High contamination</p>	 <p>1</p> <ul style="list-style-type: none"> Define a surface area of 100cm² and swab it several times with a swab. Keep the swab in a transport medium until the test performance. 															
<p>Identification</p>  <p>Salmonella spp. • positive agglutination test</p> <p>Citrobacter spp. • negative agglutination test</p> <p>Proteus spp. Providencia spp.</p> <p>Escherichia coli • positive indole test</p> <p>Pseudomonas spp. • positive oxidase test</p> <p>Listeria spp. • culture medium confirmation</p> <p>Staphylococcus aureus • positive agglutination test</p> <p>Yeasts and Moulds • confirm at the microscope</p>	 <p>2</p> <ul style="list-style-type: none"> After collection, immerse the swab in the vial of Inoculum Solution and wait 5 minutes. Distribute 0.2 mL of the sample suspension into each well of the system. 															
	 <p>3</p> <ul style="list-style-type: none"> Cover the well 4-LDC and 5-H-S with one drop of vaseline oil and incubate at 36±2°C for 18-24 hours. 															
	 <p>4</p> <ul style="list-style-type: none"> Add 2 drops of Kovacs Reagent to the well 8-IND for confirming the relevant ID. 															
	 <p>5</p> <ul style="list-style-type: none"> Confirm Yeasts and Moulds by taking a drop from the well 12-Y/M and watch for spores at the microscope (400x). 															
	 <p>6</p> <ul style="list-style-type: none"> Note down the test results on the form provided in the kit, or read and interpret with automatic readers. 															
<table border="1"> <thead> <tr> <th>DESCRIPTION</th> <th>PACKAGING</th> <th>REF.</th> </tr> </thead> <tbody> <tr> <td>HACCP SYSTEM Plus</td> <td>40 tests</td> <td>71655</td> </tr> <tr> <td>KOVACS - REAGENT</td> <td>4x25 mL</td> <td>86271</td> </tr> <tr> <td>OXIDASE TEST STICK</td> <td>50 sticks</td> <td>88029</td> </tr> <tr> <td>VASELINE OIL</td> <td>4x50 mL</td> <td>86179</td> </tr> </tbody> </table>	DESCRIPTION	PACKAGING	REF.	HACCP SYSTEM Plus	40 tests	71655	KOVACS - REAGENT	4x25 mL	86271	OXIDASE TEST STICK	50 sticks	88029	VASELINE OIL	4x50 mL	86179	
DESCRIPTION	PACKAGING	REF.														
HACCP SYSTEM Plus	40 tests	71655														
KOVACS - REAGENT	4x25 mL	86271														
OXIDASE TEST STICK	50 sticks	88029														
VASELINE OIL	4x50 mL	86179														

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

BACTERIAL COUNT

- Low microbial contamination (4-50 CFU/cm)
- Moderate microbial contamination (50-500 CFU/cm)
- High microbial contamination (>500 CFU/cm)

MICROBIAL IDENTIFICATION

- Salmonella spp.
- Citrobacter spp.
- Proteus spp./ Providencia
- Escherichia coli
- Pseudomonas spp.
- Listeria spp.
- Staphylococcus aureus
- Yeasts and Moulds

BENEFITS

- Direct inoculation from surface samples
- Conveniently designed for two tests
- Results ready in 18-24 hours
- Clear and well defined color changes
- Easy Interpretation

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
RESULTADOS

IDENTIFICADOR

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- Homogenize a quantity of foodstuff (5, 30 or 25 g) in an appropriate volume (e.g. 90 or 225 mL) of BUFFERED PEPTONE WATER or RINGER'S Solution. The final dilution of the sample should be 1:10.
- Transfer 10 mL of homogenized sample into a suitable tube. (dili. 1:10)
- Incubate at 36±1°C for 4-6 hours.

- Transfer 0.5 mL of sample suspension into a vial of Physiological Solution contained in the kit.
- Distribute 0.2 mL of suspension into each well of the system.
- Cover the wells 1-LDC, 2-H and 3-UR with one drop of vaseline oil and incubate at 36 ±1°C for 18-24 hours.

- Add 2 drops of Kovacs' Reagent to the well 8-IND for confirming the relevant ID.
- Add 2 drops of reagent H +O₂ to the well 11-CAT and watch for the formation of bubbles.
- Confirm Yeasts and Moulds by taking a drop from the well 12-Y/M and watch for spores at the microscope (400x).
- Note down the test results on the form provided in the kit, or read and interpret with automatic readers.

DESCRIPTION	PACKAGING	REF.
FOOD SYSTEM	40 tests	73690
H ₂ O ₂ REAGENT	10 mL	80027
KOVACS - REAGENT	4x25 mL	80221
OXIDASE TEST STICK	50 sticks	88019
VASELINE OIL	4x50 mL	80229

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
VALIDATED TO ISO 16140

Microbiology of food and animal feeding stuffs: protocol for the validation of alternative methods

- Listeria spp.
- Salmonella spp.

BENEFITS

- Inoculation from foodstuffs suspension
- Conveniently designed for two tests
- Results ready in 18-24 hours
- Clear and well defined color changes
- Easy interpretation

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